

## The Anfatec Level AFM – a short description

### Atomic Force Microscopy - approved devices for affordable prices

Our system is complete for almost all typical applications. It provides all basic modes as:

- ◆ high resolution contact mode
- ◆ high resolution dynamic mode (non-contact mode) with simultaneous phase acquisition
- ◆ lateral force mode (friction mode)

but also some more enhanced experiments like

- ◆ force spectroscopy
- ◆ Magnetic Force Microscopy (2<sup>nd</sup> trace imaging) – without special cantilevers & samples
- ◆ Electrical Force Microscopy (2<sup>nd</sup> trace imaging) – without special cantilevers & samples
- ◆ Nano-Lithography with script-language
- ◆ Elastic Force Microscopy (Force Modulation Mode)

and can be ordered with options for EFM (2<sup>nd</sup> 6-phase lock-in amplifier) and conduction AFM.

The instrument is a complete working AFM with full support and service!

It is as easy to handle, as a laser deflection AFM can be. We provide a fast, cost effective and uncomplicated support for our systems. With all systems, we give a personal 8 h introduction. The operation modes are demonstrated at your samples to make the start-up easy for you.

### SYSTEM PARAMETERS:

lateral resolution:	< 1 nm (practical resolution) technical resolution: 0.19 nm (18 bit achieved technical resolution) mathematical resolution: 32 Bit (< 0.1 pm)
height resolution:	< 150 pm noise floor in DNC (atomic steps and layers) technical resolution: 0.026 nm (18 bit achieved technical resolution)
maximum scan range:	50 µm (standard, others possible on request), z-range: 6 µm
maximum sample size:	4 cm x 6 cm
manual positioning range:	5 mm x 5 mm
accessories:	15 cantilevers; 1 calibration grating UMG01 20 sample holders; 2 sample boxes, tweezers

### THE SYSTEM CONSISTS OF:

- ◆ a base plate made from stone with wiring
- ◆ vibration isolation
- ◆ microscope base:
  - ☑ 3 integrated miniaturized stepper motors for head levelling
  - ☑ lateral coarse positioning with 6 mm travel range
  - ☑ self-adjustable grooves (head always in same position)
  - ☑ calibrated scanner (about 30 µm range)
  - ☑ electrical contact to sample (can be used for Electrical Force Microscopy)
- ◆ standard AFM-head
  - ☑ laser diode maximum 3 mW, 670 nm with lens system
  - ☑ laser adjustment in three axis
  - ☑ integrated 4-quadrant photo-detector with amplifier electronics



- ☑ adjustment of the laser beam onto the photo detector in two directions (X, Y)
- ☑ built-in dither piezo for acoustical excitation in dynamic mode
- ☑ integrated illumination
- ☑ color-camera with microscope optic with a direct view onto the cantilever
- ☑ uncomplicated mounting of the cantilever chips
- ◆ control electronics consisting of:
  - ☑ low-noise high voltage amplifier V45E (1 pm noise floor)
  - ☑ dual DSP controlling system incl. DS4L-Module with Interface to the **AMU 2.x**
  - ☑ control of the level station
  - ☑ 8 x 24-bit D/A and 8 x 24-bit A/D channels
  - ☑ switchable laser power
- ◆ control computer:
  - ☑ typically: up-to-date AMD-processor based modern system
  - ☑ TFT monitor(s)
  - ☑ USB video camera
  - ☑ installed software: Windows XP Prof.,  
Anfatec Scan with GNU GPL, Anfatec Present, Curve, Acrobat, OpenOffice
- ◆ Anfatec Measurement Interface **AMU 2.x** (PCI bus board with integrated LockIn amplifier)
- ◆ tools: 1 calibration grating UMG01, 1 start-up set of cantilevers (15 pcs.), connection cables
- ◆ English or German manual, certain tutorials for NC AFM, contact AFM, Force Spectroscopy, LFM, Scanner calibration
- ◆ Special features:
  - software- or hardware-linearisation and calibration of the scanner
  - mixed frequency board AMU2.x with better signal to noise ratio
  - 2<sup>nd</sup> trace imaging for MFM
  - 2<sup>nd</sup> internal LockIn amplifier for integrated Kelvin feedback and EFM imaging

## MICROSCOPE MECHANICS

The mechanics consists of three main parts: the base plate, the “body” and the “head”.

The head holds the cantilever simply by a spring loaded mechanism. It needs no glue or cantilever holder and can work with all commercially available cantilevers, even high-frequency cantilevers. All electronic components for laser and photo diode, a specialized lens and mirror system and fine mechanics are integrated in the head.

Additionally, a small CCD camera with top-view onto the tip is mounted on the head. The camera image shows directly the cantilever from the top. LED light, whose intensity is adjustable, illuminates the tip and the sample.


The base plate provides the electrical connections and allows to bring the whole system under a vacuum bell jar. This can minimize acoustical coupling from the environment and enables to work under protection gas.

The body is heavy and hangs vibration damped above the base plate. It includes all tools for coarse positioning and the scanner in 5 mm by 5 mm range. Usually, the x and y coarse movement is done manually with a travel range of 5 mm. Three stepper motors allow to position the head in three degrees of freedom.


## SAMPLE HOLDER AND SCANNER

The scanner is a tripod-type and scans the sample, while the tip is mounted in the head. The sample is mounted on a steel plate on three small magnets. Possible samples sizes are up to 2 cm by 2 cm, however, bigger samples are possible, too. There are two new linearisation modes for the scanners available: a software-mode and a hardware mode. The maximum scan range is 50  $\mu\text{m}$ .

## APPLICATIONS



← 30 nm gold clusters on a gold surface  
Scan range: 2  $\mu\text{m}$ . Height scale: 36 nm.  
*Sample: University of Karlsruhe, Dr. U. Plutowski*



atomic steps on a hydrogen passivated Si(111)-surface →  
Image range: 1  $\mu\text{m}$  x 1  $\mu\text{m}$ ; height scale: 2.6 nm  
*Sample: Hahn-Meitner Insitute Berlin, Dr. H. Angermann*

## CONTROL SYSTEM

- ◆ up-to-date PC with Windows XP UK English Professional Edition
- ◆ high-quality and high-speed LockIn-amplifier on PCI-board
- ◆ dual DSP control with a fast data interface to the SPM-control unit

## ELECTRONIC CONTROL

The system is provided with a HV amplifier [V45E](#) and a high performance SPM control unit:

- ◆ 8 channel – 24 Bit A/D (4 used for AFM, 2 additional used for hardware linearisation option)
- ◆ 8 channel – 24 Bit D/A (4 pre-specified for AFM)
- ◆ motor control
- ◆ head control
- ◆ fast digital interface to the PC


# SOFTWARE

The Software works under NT based Windows systems (typically: multilingual Windows XP).

- ◆ Image Acquisition Software "Anfatec Scan" (GPL licence)
- ◆ Image Evaluation Software "Anfatec Present"

## Anfatec Scan



- ☑ Simultaneous acquisition of up to 8 external channels + height image + LockIn channels + 8 external channels from a second basis module
- ☑ forward-trace images, backward trace images and 2<sup>nd</sup> trace images
- ☑ software feedback (PI type) with 6 different feedback modes
- ☑ free parameter input during scan
- ☑ images saved in Anfatec file format (read by Anfatec Present and SPIP<sup>(1)</sup>) and unscaled as Bitmap
- ☑ distance spectroscopy / voltage spectroscopy
- ☑ coarse positioning
- ☑ automatic approach
- ☑ switching between dynamic mode and contact mode
- ☑ adaptable to almost every hardware due to
  - ☑ free scalability of all channels in physical units
  - ☑ software offset correction
  - ☑ invert-channel-option
- ☑ user settings are saved in an initialization file



All functions are easy to select by menu buttons. A time-scaled oscilloscope provides observation of all input channels. A *Crosshairs* window (see figure) makes the adjustment of the laser for AFM easy.

## Anfatec Present

- diverse filter functions (high pass, low pass, Fourier, noise, ..)
- enhanced line and plane correction with various selection options
- data import and export
- shows scan information
- histogram for brightness & contrast
- variable colour palettes
- **3D view**



- line and **roughness analysis**
- automated unit cell detection and averaging
- copy and save functions to other programs

# LOCK-IN-AMPLIFIER ON THE AMU2.X

## Signal Input

Voltage Input	SMB
Input Impedance	1 M $\Omega$
Damage Threshold	> +/- 12 V
Bandwidth	dc to > 1 MHz (3 dB bandwidth)
Full Scale Sensitivity	7 V, 0.7 V, 70 mV
Input Noise:	
@ 100 kHz, high dynamic	< 2 $\mu V / \sqrt{Hz}$
@ 100 kHz, normal dynamic	< 0,4 $\mu V / \sqrt{Hz}$
@ 100 kHz, low dynamic	< 10 nV / $\sqrt{Hz}$

## Reference Output

Internal Oscillator	3 mHz .. 1 MHz
Frequency Resolution	3 mHz
Frequency Accuracy	+/- 50 ppm from 0°C to 70 °C
Reference Output Voltage	< 1 mVpp ... max. 20 Vpp
Output Noise @ 100 kHz for 7 V <sub>rms</sub> output	160 nV / $\sqrt{Hz}$

## SCANNER:

Maximum range in z-direction:	6 $\mu m$ +/- 0.3 $\mu m$ / 150V
Maximum range in x- and y-direction:	50 $\mu m$
technical resolution in x- und y-direction:	0,9 nm
technical resolution in z-direction:	0,034 nm
achieved resolution:	about 5 nm
achieved resolution z-direction:	< 0.2 nm

## OPTIONAL FEATURES (NOT NECESSARY FOR STANDARD APPLICATIONS):

- ◆ Vibration isolation table under the microscope
- ◆ Hardware scanner linearisation
- ◆ Glass bell jar for acoustic protection
- ◆ Additional cantilever packages and gratings
- ◆ enhanced LFM mode sensitivity due to a spot-like laser diode
- ◆ additional LockIn amplifier for dynamic EFM or MFM
- ◆ implemented Kelvin feedback
- ◆ current amplifier for conductance AFM incl. power supply
- ◆ SPIP – Scanning Probe Image Processor - with all customer specific modules from Imaging Metrology
- ◆ 2<sup>nd</sup> TFT monitor